FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Material Safety Data Sheet (MSDS)

MSDS #: 571.10

Revision Date: November 25, 2002

Section 1 — Chemical Product and Company Identification

Oxalic Acid Solution

Flinn Scientific, Inc. P.O. Box 219 Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300

Section 2 — Composition, Information on Ingredients

Oxalic Acid (6153-56-6) 40%, and Water (7732-18-5) 60%.

CAS#: None established

Section 3 — Hazards Identification

Clear solution with white solid sitting on the bottom of the container. Odorless. Moderately toxic by ingestion. Corrosive to body tissues. Avoid all body tissue contact.

FLINN AT-A-GLANCE
Health-2
Flammability-0
Reactivity-1
Exposure-2

Storage-0

0 is low hazard, 3 is high hazard

Section 4 — First Aid Measures

Call a physician, seek medical attention for further treatment, observation and support after first aid.

Inhalation: Remove to fresh air at once. If breathing has stopped give artificial respiration immediately.

Eye: Immediately flush with fresh water for 15 minutes.

External: Wash continuously with fresh water for 15 minutes.

Internal: Give 1 to 2 cups of water or milk, followed by a gastric antacid, such as milk of magnesia. Do not induce vomiting. Call a physician or poison control at once.

Section 5 — Fire Fighting Measures

Non flammable, non combustible liquid.

Decomposition releases toxic fumes.

NFPA CODE
None Established

Fire Fighting Instructions: Use triclass, dry chemical fire extinguisher. Firefighters should wear PPE and SCBA with full facepiece operated in positive pressure mode.

and SCBA with full racepiece operated in positive pressure mode.

Section 6 — Accidental Release Measures

Restrict unprotected personnel from area and ventilate area. Contain spill with sand or absorbent material; deposit in sealed bag or container. See Sections 8 and 13 for further information.

Section 7 — Handling and Storage

Flinn Suggested Chemical Storage Pattern: Organic #1. Store with acids, anhydrides and peracids.

<u>Section 8 — Exposure Controls , Personal Protection</u>

Avoid contact with eyes, skin and clothing. Wear chemical splash goggles, chemical-resistant gloves and chemical-resistant apron.

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Material Safety Data Sheet (MSDS)

Oxalic Acid Solution

MSDS #: 571.10

Revision Date: November 25, 2002

Section 9 — Physical and Chemical Properties

Clear solution with white solid sitting on the bottom of the container.

Section 10 — Stability and Reactivity

Avoid contact with bases, acid chlorides, alkali metals. Shelf life: Indefinite.

Section 11 — Toxicological Information

Acute effects: Irritant, toxic Chronic effects: N.A. Target organs: N.A. ORL-RAT LD50: 375 mg/kg as oxalic acid.

IHL-RAT LC50: N.A. SKN-RBT LD50: N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

Section 12 — Ecological Information

Data not yet available.

Section 13 — Disposal Considerations

Please consult with state and local regulations. Flinn Suggested Disposal Method #24a is one option.

Section 14 — Transport Information

Shipping Name: Not regulated

Hazard Class: N/A UN Number: N/A N/A = Not applicable

Section 15 — Regulatory Information

Not listed.

Section 16 — Other Information

Consult your copy of the Flinn Scientific Catalog/Reference Manual for additional information about laboratory chemicals. This Material Safety Data Sheet (MSDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. Flinn Scientific Inc. assumes no legal responsibility for use or reliance upon this data.

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Questions on Chemical Disposal or Storage?--Call Flinn

flinn@flinnsci.com www.flinnsci.com P.O. Box 219 Batavia IL 60510 (800) 452-1261 Fax (866) 452-1436